

DELIVERY SYSTEMS & LIGHTING

50
YEARS

Proma has an extensive history of dental delivery innovation. Proma first introduced the self-contained patient water system, ambidextrous delivery of doctor and assistant's instrumentation, and pinch technology to control handpiece air and water.



Proma

betterDentistry | byDesign

www.royaldentalgroup.com

ROYAL
DENTAL GROUP

MOUNTING STYLES

Compass

The Compass delivery is Prisma's most popular system, providing over-the-patient positioning of instruments with left/right flexibility for multiple operators.


- A6210 **DOCTOR** (shown)
- A6220 **EURO**
- A6230 **HYGIENE**
- A6380 **DUAL - VALUE LINE**



North

The Prisma North delivery provides 180° arrangement of assistant's or hygiene instrumentation, mounted at the 12 o'clock position on any Royal chair.

- A6410 **HYGIENE**
- A6760 **ASSISTANT** (shown)



Mobile

Proma's mobile delivery systems move freely throughout the operatory for ultimate flexibility in workflow and maximum utilization of space.

A6612 **DOCTOR** (shown)

A6622 **HYGIENE**

A6632 **DUAL**

A6642 **ASSISTANT**



Traditional Post


The post mounted, over-the-patient delivery provides flexibility while utilizing minimal operatory space. Proma post systems can be installed on nearly all dental chairs.

A6310 **DOCTOR** (shown)

A6721 **ASSISTANT w/ CUSPIDOR** (shown)

A6330 **DOCTOR - VALUE LINE**

A6340 **DUAL - VALUE LINE**



Floor

Floor mounted systems easily integrate with custom cabinetry and uniquely designed treatment rooms.

A6552 **DUAL** (shown)

A6554 **HYGIENE**

A5562 **ASSISTANT**



DELIVERY CONFIGURATIONS



DOCTOR
3 handpieces
Air/water syringe



DOCTOR EURO
3 handpieces
Air/water syringe



ASSISTANT
Air/water syringe
HVE
Saliva ejector



HYGIENE
2 handpieces
Air/water syringe
HVE
Saliva ejector



DUAL DR + ASST
3 handpieces
2 Air/water syringes
HVE
Saliva ejector

LIGHTING



REFLECTED LED

Proma's LED light utilizes dual parabolic lenses for reflective illumination of the oral cavity. The color corrected pattern is brightest in the center, soft around the edges and nearly shadow-free.



HALOGEN

Proma's traditional halogen light uses a quartz halogen bulb for color corrected illumination, which has been the dependable mainstay of dentistry lighting for decades.

MOUNTING OPTIONS:



A5110LED **CHAIR**



A5150LED **TRACK**



A5120LED **CEILING**



A5140LED **CABINET**

DEPENDABILITY

The delivery system is the financial foundation of your practice. Flexible positioning that adapts to each member of your staff maximizes productivity and facilitates healthy ergonomic posture. Rely on Proma units to provide years of reliable service.



Proma